Shake no more All-in-one

free from manual shaking

Free from Physical burden



Set then press a button

No need to shake a container for a given amount of time with force.

Safe measuring environment

No accidents caused by the unit slipping from your hands while shaking.

02 Free from Inconsistent measurement value

03 Free from





Consistent measurement results regardless of measurement taker

Compared to the "manual shaking" method, measurement results from an automatic measurement will not be affected by the measurement taker.

No cumbersome calculation after measurement


Using the preexisting and current data obtained by CooRe, a calibration curve can be created which can then be used to program a coefficient. Measurement values are automatically computed based off the past data, which means no stored data is wasted.

~Company "S" Quality Assurance Department~ Unlike traditional instruments, it does not require intensive swinging motions, instead it can be readily used.

~Company "I" Product Development Department~

There is no more variation in measurement values that depend on who is measuring the sample.

Name of Compartment



Cat.No.	9332	9334	Temperature	Brix	: 5.0 to 30.0	
Model	CooRe	CooRe 3L	Compensation Range	CO ₂ Concen	tration: 0 to 25°C	: 0 to 25 °C
Measurement Scales	Brix (Automatic temperature compensation) CO2 Concentration (Automatic temperature compensation) Pressure Temperature		Minimum Sample Volume	Approx. 100 mL *1st measurement 200mL (Pressure chamber volume : 20mL)		
			Data Display History	Number Measurement year month date time		
Measurement System	CO ₂ Concentration	: Refractive index method : Converted from pressure : Sample agitation measurement method		Scale Brix Temperature when Brix is measured COs concentration		
Measurement Range		: 0.00 to 20.00 % : 0.000 to 12.000 vol. 0.000 to 24.000 × 10 ⁻⁶ kg/cm ³		Temperature when CO ₂ concentration is measured Pressure Temperature when pressure is measured		
		0.000 to 24.000 g/L : 0.00 to 10.00 bar 10 to 10000 mbar 0.000 to 145.000 psi 0.000 to 1.000 MPa 0 to 1000 KPa : 0.0 to 30.0 °C 32.0 to 86.0 °F 273 to 303 K	Output	USB memory (CSV text file format)		
			Pressure Resistance	Pressure cha	mber : 1Mpa : 1MPa	
			International Protection Class	IP65 *Applicable only if the USB terminal cover is on a power terminal cover is on (or AC power is control or AC p		
			Power Supply	AC 100 to 240 V ,50/60Hz		
			Power Consumption	70 VA		
Resolution	Brix : 0.01 % CO ₂ Concentration : 0.001 vol. 0.001 × 10 ⁶ kg/cm ³ 0.001 g/L		Dimensions & Weight	CooRe 3L	Approx. 40.5 × 21: 29lbs/13kg (Power Approx. 40.5 × 21: 29lbs/13kg (Power	cable included.) × 51.5 cm
	Pressure	: 0.01 bar 10 mbar 0.001 psi 0.001 MPa 1 kPa	Content			

: 0.1 °C 0.1 °F

: ±1.0 % (at F.S. 20 °C)

Brix (Absolute value) : ±0.05 %

Output	USB memory (CSV text file format)		
Pressure Resistance	Pressure ch Prism	namber : 1Mpa : 1MPa	
International Protection Class		icable only if the USB terminal cover is on and er terminal cover is on (or AC power is connected).	
Power Supply	AC 100 to 240 V ,50/60Hz		
Power Consumption	70 VA		
Dimensions & Weight	CooRe	Approx. 40.5 × 21 × 45 cm 29lbs/13kg (Power cable included.)	
	CooRe 3L	Approx. $40.5 \times 21 \times 51.5$ cm 29 lbs/13kg (Power cable included.)	
Content			
1. Main unit·····			


Part Name	Part No.		
Lithium-ion Battery	RE-74840		
Pulp Filter Unit	RE-79424		
Filter Cartridge (for Pulp Filter Unit)	RE-79425		
Desiccant (100g)	RE-99440		
Screw set for Back Cover	RE-99441		
Extension Drain Tube	RE-78068		
Shaft Lubricant for CooRe	RE-99442		

3. Extension Drain Tube ··

6. AC power cable...

5. Shaft Lubricant for CooRe

4. Connector (for Extension Drain Tube)



All ATAGO products are designed and manufactured in Japan.



TEL: 81-3-3431-1943 FAX: 81-3-3431-1945



ATAGO products comply with HACCP,GMP,



Djupdalsvägen 18 S-192 51 SOLLENTUNA Telefon: 08-83 08 80 Telefax: 08-82 33 44 E-post: info@bohlins-maskiner.se


* Specifications and appearance are subject to change without notice.

Copyright © 2019 ATAGO CO., LTD. All rights reserved. ENV.04 1906800PP Printed in Japan
















Readily take measurements any time, anywhere

No other automatic gas volume meter allows taking measurements any time and anywhere. CooRe is equipped with many convenient functions and components.

World first all-in-one piercing and measuring instrument

The piercing and measuring sections are together in one unit. Traditional automatic gas volume meters are equipped with nitrogen gas cylinders which weigh almost 10kg/22lbs. With CooRe, the pump is internalized requiring no need to connect or change the nitrogen gas cylinder.

Excellent operability

Easy piercing

This instrument can pierce any type of bottle easily and smoothly, whether can or bottle. In addition, it is equipped with a cover to prevent liquid from scattering when a bottle is removed from the instrument.

Sample amount is 100mL

Only 100mL of sample is required for measurements. If the capacity is that of the common soda bottle, all containers can be measured.

Cover

Cover to prevent liquid or pieces of a damaged can, glass bottle, or PET bottle from scattering.

Portable

Can be carried with one hand

The unit approximately weighs 13kg/29lbs. The weight being about half of traditional automatic gas volume meters, it is very easy to handle. The display and measurement sections are all in one place with a handle on the top to free up one of your hands during transit.



Can be used in area with no power supply

The power source can either be AC adapter or lithium ion battery. No more thinking about where to measure. Moreover, the rechargeable lithium ion battery can be used repeatedly, reducing running costs.

One minute to measure

After pressing the button, the measurement will automatically display in one minute. *1st measurement approx. 110 seconds

ATAGO's Reliability


ATAGO is a Brix meter manufacturer since 1940 serving 154 countries worldwide.






Touchscreen

Smooth and absolute operation



Measurement Screen

Before taking a measurement, select "LEAK TEST" then press "OK" once the completion message is displayed. Next, press "MEASURE" to begin taking measurements.



Value Correction Input Screen

Offset feature


Using the preexisting and current data obtained by CooRe, a calibration curve can be created which can then be used to program a coefficient.

Measurement values are automatically computed based off of past data, which means no stored data is wa




Scale List Screen

By programming the scale, there is no need to reset the unit or mode



Use this screen



™ Data History



Settings Screen



Gateway to access basic configuration Navigate to each configuration page using the icons. Confusion-free, intuitive operation.







